Cistatin C

Uputa o lijeku za imunoanalizu cistatina C na sustavima serije AU društva Beckman Coulter® (AU5800, AU680, AU480, AU2700, DxC 700 AU)

REF B08179

Namjena

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Reagens sustava: B08179

ID reagensa: 228

Specific Test Parameters										
General		LIH	ISE	Range						
Test Name:		CysC	<	>	Type:	Serum	Operation:	Yes		
Sample:	Volume	2	µL	Dilution	0	µL	Pre-Dilution Rate:	1		
Reagents:	R1 Volume	150	µL	Dilution	0	µL	Min OD		Max OD	
	R2 Volume	30	µL	Dilution	10	µL	L		H	
Wavelength:	Pri.	540	▽	Sec.		▽	Reagent OD limit:			
Method:	End Point		▽				First L	-2.0	First H	2.0
Reaction slope:		+	▽				Last L		Last H	
Measuring Point 1:	First	13		Last	27		Dynamic Range:			
Measuring Point 2:	First			Last			L	0.4	H	7.8
Linearity:			%				Correlation Factor:			
No Lag Time:			▽				A	1.00	B	0.00
							On-board stability period:	60		

Specific Test Parameters									
General		LIH	ISE	Range					
Test Name:		CysC	<	>	Type:	Serum			
Value/Flag:	#	▽	Level L:	#	Level H:	#			
Normal Ranges:		Age L	Age H						
	Sex	Year	Month	Year	Month	L	H		
<input type="checkbox"/>	1.	#	▽	#	#	#	#		
<input type="checkbox"/>	2.	#	▽	#	#	#	#		
<input type="checkbox"/>	3.	#	▽	#	#	#	#		
<input type="checkbox"/>	4.	#	▽	#	#	#	#		
<input type="checkbox"/>	5.	#	▽	#	#	#	#		
<input type="checkbox"/>	6.	#	▽	#	#	#	#		
	7.	None Selected				#	#		
	8.	Out of Range		L	H	#	#		
Panic Value:	#		#		#	Unit:	mg/L	Decimal places:	#

Calibration Specific											
General		ISE									
Test Name:		CysC	<	>			Type:	Serum			
Calibration Type:		6AB	▽	Formula:	Spline	▽	Counts:	#	Process:	CONC	▽
Point 1:	Cal. No.	1	OD	CONC	**	Factor/OD-L		Factor/OD-H			
Point 2:	2			**							
Point 3:	3			**							
Point 4:	4			**							
Point 5:	5			**							
Point 6:	6			**							
Point 7:				**							
1-Point Cal. Point:	<input type="checkbox"/>	with CONC-0	Slope Check:	+	▽	Advanced Calibration:	#	▽			
MB Type Factor:		Calibration Stability Period:	28								

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Reagens sustava: B08179

ID reagensa: 228

Parameters		Specific Test Parameters			
General	LIH	ISE	Calculated Test	Range	
Test Name: <input type="text" value="CysC"/> < >		Type: <input type="text" value="Serum"/>		Operation: <input type="text" value="Yes"/>	
Sample Volume	<input type="text" value="2"/> μL	Dilution	<input type="text" value="0"/> μL	OD Limit	
Pre-Dilution Rate	<input type="text" value="1"/>	Min.OD		Max.OD	
Rgt. Volume R1(R1-1)	<input type="text" value="150"/> μL	Dilution	<input type="text" value="0"/> μL	Reagent OD Limit	
R2(R2-1) <input type="text" value="30"/> μL		Dilution	<input type="text" value="10"/> μL	First	High
Common Rgt. Type		Name		Low	<input type="text" value="-2.0"/>
Wavelength	Pri <input type="text" value="540"/> nm	Sec. <input type="text" value=""/>	Dynamic Range Low	<input type="text" value="0.4"/>	High <input type="text" value="7.8"/>
Method	<input type="text" value="End Point"/>	Correlation Factor A		<input type="text" value="1.00"/>	B <input type="text" value="0.00"/>
Reaction Slope	<input type="text" value=""/>	Factor for Maker A		<input type="text" value=""/>	B <input type="text" value=""/>
Measuring Point1 First	<input type="text" value="13"/>	Last	<input type="text" value="27"/>	Onboard Stability Period <input type="text" value="60"/> Day <input type="text" value=""/> Hour	
Measuring Point2 First	<input type="text" value=""/>	Last	<input type="text" value=""/>	LIH Influence Check <input type="text" value=""/>	
Linearity Limit	<input type="text" value=""/>	Lipemia		<input type="text" value=""/>	
Lag Time Check	<input type="text" value=""/>	Icterus		<input type="text" value=""/>	
		Hemolysis		<input type="text" value=""/>	

Parameters		Specific Test Parameters					
General	LIH	ISE	Calculated Test	Range			
Test Name: <input type="text" value="CysC"/> < >		Type: <input type="text" value="Serum"/>					
Value/Flag:	<input type="text" value="#"/>	Low	<input type="text" value=""/>	High	<input type="text" value=""/>		
Level	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>		
Specific Ranges:	From	To	Low	High	Panic Value		
	Sex	Year	Month	Year	Month	Low	High
<input type="checkbox"/> 1.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
<input type="checkbox"/> 2.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
<input type="checkbox"/> 3.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
<input type="checkbox"/> 4.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
<input type="checkbox"/> 5.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
<input type="checkbox"/> 6.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
7. No demographics						<input type="text" value=""/>	<input type="text" value=""/>
8. Not within expected values						<input type="text" value=""/>	<input type="text" value=""/>
Unit		<input type="text" value="mg/L"/>	Decimal Places		<input type="text" value=""/>		

Parameters		Calibration Parameters				
Calibrators	Calibration Specific	STAT Table Calibration				
General	ISE					
Test Name: <input type="text" value="CysC"/> < >		Type: <input type="text" value="Serum"/>		<input type="checkbox"/> Use Serum Cal.		
Calibration Type: <input type="text" value="6AB"/>		Formula: <input type="text" value="Spline"/>		Counts: <input type="text" value=""/>		
<Calibrator Parameters>		Factor Range		Slope Check <input type="text" value="+"/> <		
Calibrator	OD	Conc	Low	High	Allowance Range Check	
Point 1:	<input type="text" value="1"/>	<input type="text" value="**"/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="checkbox"/> Reagent Blank <input type="text" value=""/>	
Point 2:	<input type="text" value="2"/>	<input type="text" value="**"/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="checkbox"/> Calibration <input type="text" value=""/>	
Point 3:	<input type="text" value="3"/>	<input type="text" value="**"/>	<input type="text" value=""/>	<input type="text" value=""/>	Advanced Calibration	
Point 4:	<input type="text" value="4"/>	<input type="text" value="**"/>	<input type="text" value=""/>	<input type="text" value=""/>	Operation <input type="text" value=""/> <	
Point 5:	<input type="text" value="5"/>	<input type="text" value="**"/>	<input type="text" value=""/>	<input type="text" value=""/>	Interval (RB/ACAL) <input type="text" value=""/> <	
Point 6:	<input type="text" value="6"/>	<input type="text" value="**"/>	<input type="text" value=""/>	<input type="text" value=""/>	<Point Cal. For No. of Correction Points <input type="text" value=""/>	
Point 7:	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	Use Master Curve <input type="checkbox"/> <	
Point 8:	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="checkbox"/> Lot Calibration	
Point 9:	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	Master Curve>	
Point 10:	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	Calibrator	
	OD	Conc	Low	High	Stability	
Point-1	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	Reagent Blank <input type="text" value="28"/> Day <input type="text" value=""/> Hour	
Point-2	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	Calibration <input type="text" value="28"/> Day <input type="text" value=""/> Hour	
MB Type Factor: <input type="text" value=""/>		1-Point Calibration Point <input type="text" value=""/>		<input type="checkbox"/> with Conc-0		

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Reagens sustava: B08179
Naziv pretrage: CYS1G

ID reagensa: 228

General	LIH	ISE	Calculated Test	Range
Test Name: CYS1G <input type="button" value="Test No"/> Type: Serum <input type="button" value="Operation"/> Yes				
Sample Volume	<input type="text" value="2.0"/> µL	Dilution	<input type="text" value="0"/> µL	OD Limit
Pre-Dilution Rate	<input type="text" value="1"/>			Min. OD <input type="text"/> Max OD <input type="text"/>
Reagent Volume	R1 (R1-1) <input type="text" value="150"/> µL	Dilution	<input type="text" value="0"/> µL	Reagent OD Limit 1 st Low <input type="text" value="-2.0000"/> High <input type="text" value="2.0000"/>
	R1-2 <input type="text"/> µL	Dilution	<input type="text"/> µL	Last Low <input type="text" value="-2.0000"/> High <input type="text" value="2.0000"/>
	R2 (R2-1) <input type="text" value="30"/> µL	Dilution	<input type="text" value="10"/> µL	Analytical Measuring Range Low <input type="text" value="0.4"/> High <input type="text" value="7.8"/>
Common Reagent	Type <input type="text" value="None"/>	Name	<input type="text" value="None"/>	Correlation Factor A <input type="text" value="1"/> B <input type="text" value="0"/>
Wavelength	Pri <input type="text" value="540"/> nm	Sec	<input type="text" value="None"/> nm	Manufacturer Factor A <input type="text" value="1"/> B <input type="text" value="0"/>
Method	<input type="text" value="END"/>			
Reaction Slope	<input type="text" value="+"/> ▼			Onboard Stability Period <input type="text" value="60"/> Day <input type="text" value="0"/> Hour
Measuring Point-1	1st <input type="text" value="13"/>	Last	<input type="text" value="27"/>	LIH Influence Check <input type="text" value="No"/> ▼
Measuring Point-2	1st <input type="text"/>	Last	<input type="text"/>	Lipemia <input type="text"/> ▼
Linearity Limit	<input type="text"/> %			Icterus <input type="text"/> ▼
Lag Time Check	<input type="text"/>			Hemolysis <input type="text"/> ▼

General	LIH	ISE	Calculated Test	Range
Test Name: CYS1G <input type="button" value="Test No"/> Type: Serum				
Value/Flag	<input type="text" value="Value"/>	Level	Low <input type="text" value="-99999.99"/> High <input type="text" value="99999.99"/>	
Specific Ranges	From	To		
<input type="checkbox"/> 1:	Sex <input type="text" value="#"/> ▼	Year <input type="text" value="#"/> Month <input type="text" value="#"/>	Year <input type="text" value="#"/> Month <input type="text" value="#"/>	Other Type <input type="text" value="None"/> Low <input type="text" value="#"/> High <input type="text" value="#"/>
<input type="checkbox"/> 2:	<input type="text" value="#"/> ▼	<input type="text" value="#"/> <input type="text" value="#"/>	<input type="text" value="#"/> <input type="text" value="#"/>	<input type="text" value="None"/> <input type="text" value="#"/> <input type="text" value="#"/>
<input type="checkbox"/> 3:	<input type="text" value="#"/> ▼	<input type="text" value="#"/> <input type="text" value="#"/>	<input type="text" value="#"/> <input type="text" value="#"/>	<input type="text" value="None"/> <input type="text" value="#"/> <input type="text" value="#"/>
<input type="checkbox"/> 4:	<input type="text" value="#"/> ▼	<input type="text" value="#"/> <input type="text" value="#"/>	<input type="text" value="#"/> <input type="text" value="#"/>	<input type="text" value="None"/> <input type="text" value="#"/> <input type="text" value="#"/>
<input type="checkbox"/> 5:	<input type="text" value="#"/> ▼	<input type="text" value="#"/> <input type="text" value="#"/>	<input type="text" value="#"/> <input type="text" value="#"/>	<input type="text" value="None"/> <input type="text" value="#"/> <input type="text" value="#"/>
<input type="checkbox"/> 6:	<input type="text" value="#"/> ▼	<input type="text" value="#"/> <input type="text" value="#"/>	<input type="text" value="#"/> <input type="text" value="#"/>	<input type="text" value="None"/> <input type="text" value="#"/> <input type="text" value="#"/>
7:	Standard demographics			<input type="text" value="#"/> <input type="text" value="#"/>
8:	Not within expected values			<input type="text" value="#"/> <input type="text" value="#"/>
Critical Limits	Low <input type="text" value="#"/>	High <input type="text" value="#"/>	Unit <input type="text" value="mg/L"/>	Select <input type="text"/> Decimal Places <input type="text" value="2"/>

Calibrators	General	ISE		
Test Name: CYS1G <input type="button" value="Use Serum Cal."/> Type: Serum				
Calibration Type:	<input type="text" value="6AB"/> ▼	Formula: <input type="text" value="Spline"/> ▼		
		Counts: <input type="text" value="2"/>		
<Calibrator Parameters>		Slope Check <input type="text" value="+"/> ▼		
	Calibrator	OD		
		Conc		
		Range		
		Low		
		High		
Point-1	CYSC Calibrator Level 1 ▼	†	-2.0000	2.0000
Point-2	CYSC Calibrator Level 2 ▼	†	-2.0000	2.0000
Point-3	CYSC Calibrator Level 3 ▼	†	-2.0000	2.0000
Point-4	CYSC Calibrator Level 4 ▼	†	-2.0000	2.0000
Point-5	CYSC Calibrator Level 5 ▼	†	-2.0000	2.0000
Point-6	CYSC Calibrator Level 6 ▼	†	-2.0000	2.0000
Point-7	<input type="text"/>			
MB Type Factor	<input type="text"/>	1-Point Calibration Point <input type="text" value="None"/> ▼	<input type="checkbox"/> with Conc-0	Interval (RB) <input type="text"/> ▼
				Interval (ACAL) <input type="text"/> ▼
			Stability	
			Reagent Blank	<input type="text" value="28"/> Day <input type="text" value="0"/> Hour
			Calibration	<input type="text" value="28"/> Day <input type="text" value="0"/> Hour

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