Cystatin C

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Rendszerreagens: B08179

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Specific Test Parameters									
General		LIH	ISE	Range					
Test Name:	CysC	<	>	Type:	Serum	Operation:	Yes		
Sample:	Volume	2	µL	Dilution	0	µL	Pre-Dilution Rate:	1	
Reagents:	R1 Volume	150	µL	Dilution	0	µL	Min OD		Max OD
	R2 Volume	30	µL	Dilution	10	µL	L		H
Wavelength:	Pri.	540	nm	Sec.			Reagent OD limit:		
Method:	End Point						First L	-2.0	First H
Reaction slope:	+						Last L		Last H
Measuring Point 1:	First	13		Last	27		Dynamic Range:		
Measuring Point 2:	First			Last			L	0.4	H
Linearity:			%				Correlation Factor:		
No Lag Time:							A	1.00	B
							On-board stability period:	60	

Specific Test Parameters									
General		LIH	ISE	Range					
Test Name:	CysC	<	>	Type:	Serum				
Value/Flag:	#	Level L:	#	Level H:	#				
Normal Ranges:	Sex	Age L	Year	Month	Age H	Year	Month	L	H
<input type="checkbox"/>	1.	#	#	#	#	#	#	#	#
<input type="checkbox"/>	2.	#	#	#	#	#	#	#	#
<input type="checkbox"/>	3.	#	#	#	#	#	#	#	#
<input type="checkbox"/>	4.	#	#	#	#	#	#	#	#
<input type="checkbox"/>	5.	#	#	#	#	#	#	#	#
<input type="checkbox"/>	6.	#	#	#	#	#	#	#	#
	7.	None Selected							
	8.	Out of Range							
			L		H				
Panic Value:	#	#	#	#	#	Unit:	mg/L	Decimal places:	#

Calibration Specific									
General		ISE							
Test Name:	CysC	<	>			Type:	Serum		
Calibration Type:	5AB	Formula:	Spline	Counts:	#	Process:	CONC		
Point 1:	Cal. No.	1	OD	CONC	**	Factor/OD-L		Factor/OD-H	
Point 2:	2			**					
Point 3:	3			**					
Point 4:	4			**					
Point 5:	5			**					
Point 6:	6			**					
Point 7:				**					
1-Point Cal. Point:	<input type="checkbox"/>	with CONC-0	Slope Check:	+	Advanced Calibration:	#			
MB Type Factor:		Calibration Stability Period:	28						

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Rendszerreagens: B08179

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Parameters		Specific Test Parameters			
General	LIH	ISE	Calculated Test	Range	
Test Name: <input type="text" value="CysC"/> < > Type: <input type="text" value="Serum"/> Operation <input type="text" value="Yes"/>					
Sample Volume	<input type="text" value="2"/> μL	Dilution	<input type="text" value="0"/> μL	OD Limit	
Pre-Dilution Rate	<input type="text" value="1"/>			Min.OD	<input type="text"/>
Rgt. Volume	R1(R1-1) <input type="text" value="150"/> μL	Dilution	<input type="text" value="0"/> μL	Reagent OD Limit	<input type="text"/>
				First	Low <input type="text" value="-2.0"/> High <input type="text" value="2.0"/>
				Last	Low <input type="text"/> High <input type="text"/>
				Dynamic Range Low	0.4 High <input type="text" value="7.8"/>
R2(R2-1)	<input type="text" value="30"/> μL	Dilution	<input type="text" value="10"/> μL	Correlation Factor A	1.00 B <input type="text" value="0.00"/>
Common Rgt. Type	<input type="text"/>	Name	<input type="text"/>	Factor for Maker A	B <input type="text"/>
Wavelength	Pri <input type="text" value="540"/> nm	Sec.	<input type="text"/>		
Method	<input type="text" value="End Point"/>				
Reaction Slope	<input type="text" value="+"/> ∇			Onboard Stability Period	<input type="text" value="60"/> Day <input type="text"/> Hour
Measuring Point1	First <input type="text" value="13"/>	Last	<input type="text" value="27"/>	LIH Influence Check	<input type="text"/>
Measuring Point2	First <input type="text"/>	Last	<input type="text"/>	Lipemia	<input type="text"/>
Linearity Limit	<input type="text"/>			Icterus	<input type="text"/>
Lag Time Check	<input type="text"/>			Hemolysis	<input type="text"/>

Parameters		Specific Test Parameters						
General	LIH	ISE	Calculated Test	Range				
Test Name: <input type="text" value="CysC"/> < > Type: <input type="text" value="Serum"/>								
Value/Flag:	<input type="text" value="#"/>	Low	<input type="text"/>	High	<input type="text"/>			
Level	<input type="text" value="#"/>							
Specific Ranges:	From	To			Panic Value			
	Sex	Year	Month	Year	Month	Low	High	
<input type="checkbox"/> 1.	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	
<input type="checkbox"/> 2.	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	
<input type="checkbox"/> 3.	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	
<input type="checkbox"/> 4.	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	
<input type="checkbox"/> 5.	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	
<input type="checkbox"/> 6.	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	
7.	No demographics						<input type="text" value="#"/>	<input type="text" value="#"/>
8.	Not within expected values						<input type="text" value="#"/>	<input type="text" value="#"/>
Unit	<input type="text" value="mg/L"/>	Decimal Places	<input type="text" value="#"/>					

Parameters		Calibration Parameters			
Calibrators	Calibration Specific	STAT Table Calibration			
General	ISE				
Test Name: <input type="text" value="CysC"/> < > Type: <input type="text" value="Serum"/> <input type="checkbox"/> Use Serum Cal.					
Calibration Type: <input type="text" value="6AB"/> Formula: <input type="text" value="Spline"/> Counts: <input type="text" value="#"/>					
<Calibrator Parameters>					
	Calibrator	OD	Conc	Factor Range	Slope Check
				Low	High
Point 1:	<input type="text" value="1"/>	<input type="text"/>	<input type="text" value="**"/>	<input type="text"/>	<input type="text"/>
Point 2:	<input type="text" value="2"/>	<input type="text"/>	<input type="text" value="**"/>	<input type="text"/>	<input type="text"/>
Point 3:	<input type="text" value="3"/>	<input type="text"/>	<input type="text" value="**"/>	<input type="text"/>	<input type="text"/>
Point 4:	<input type="text" value="4"/>	<input type="text"/>	<input type="text" value="**"/>	<input type="text"/>	<input type="text"/>
Point 5:	<input type="text" value="5"/>	<input type="text"/>	<input type="text" value="**"/>	<input type="text"/>	<input type="text"/>
Point 6:	<input type="text" value="6"/>	<input type="text"/>	<input type="text" value="**"/>	<input type="text"/>	<input type="text"/>
Point 7:	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
Point 8:	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
Point 9:	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
Point 10:	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
Allowance Range Check					
<input type="checkbox"/> Reagent Blank <input type="text"/>					
<input type="checkbox"/> Calibration <input type="text"/>					
Advanced Calibration Operation <input type="text" value="#"/>					
Interval (RB/ACAL) <input type="text" value="#"/>					
<Point Cal. For Master Curve>					
	Calibrator	OD	Conc	OD Range	Stability
				Low	High
Point-1	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	Reagent Blank <input type="text" value="28"/> Day <input type="text"/> Hour
Point-2	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	Calibration <input type="text" value="28"/> Day <input type="text"/> Hour
MB Type Factor: <input type="text"/> 1-Point Calibration Point <input type="text"/> <input type="checkbox"/> with Conc-0					

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General	LIH	ISE	Calculated Test	Range
Test Name: CYS1G <input type="button" value="Test No"/> Type: Serum <input type="button" value="Operation"/> Yes				
Sample Volume	<input type="text" value="2.0"/> µL	Dilution	<input type="text" value="0"/> µL	OD Limit
Pre-Dilution Rate	<input type="text" value="1"/>			Min. OD <input type="text"/> Max OD <input type="text"/>
Reagent Volume	R1 (R1-1) <input type="text" value="150"/> µL	Dilution	<input type="text" value="0"/> µL	Reagent OD Limit 1 st Low <input type="text" value="-2.0000"/> High <input type="text" value="2.0000"/>
	R1-2 <input type="text"/> µL	Dilution	<input type="text"/> µL	Last Low <input type="text" value="-2.0000"/> High <input type="text" value="2.0000"/>
	R2 (R2-1) <input type="text" value="30"/> µL	Dilution	<input type="text" value="10"/> µL	Analytical Measuring Range Low <input type="text" value="0.4"/> High <input type="text" value="7.8"/>
Common Reagent	Type <input type="text" value="None"/>	Name	<input type="text" value="None"/>	Correlation Factor A <input type="text" value="1"/> B <input type="text" value="0"/>
Wavelength	Pri <input type="text" value="540"/> nm	Sec	<input type="text" value="None"/> nm	Manufacturer Factor A <input type="text" value="1"/> B <input type="text" value="0"/>
Method	<input type="text" value="END"/>			
Reaction Slope	<input type="text" value="+"/>			Onboard Stability Period <input type="text" value="60"/> Day <input type="text" value="0"/> Hour
Measuring Point-1	1st <input type="text" value="13"/>	Last	<input type="text" value="27"/>	LIH Influence Check <input type="text" value="No"/>
Measuring Point-2	1st <input type="text"/>	Last	<input type="text"/>	Lipemia <input type="text"/>
Linearity Limit	<input type="text"/> %			Icterus <input type="text"/>
Lag Time Check	<input type="text"/>			Hemolysis <input type="text"/>

General	LIH	ISE	Calculated Test	Range
Test Name: CYS1G <input type="button" value="Test No"/> Type: Serum				
Value/Flag	<input type="text" value="Value"/>	Level	Low <input type="text" value="-99999.99"/> High <input type="text" value="99999.99"/>	
Specific Ranges	From	To		
<input type="checkbox"/> 1:	Sex <input type="text" value="#"/>	Year <input type="text" value="#"/> Month <input type="text" value="#"/>	Year <input type="text" value="#"/> Month <input type="text" value="#"/>	Other Type <input type="text" value="None"/> Low <input type="text" value="#"/> High <input type="text" value="#"/>
<input type="checkbox"/> 2:	<input type="text" value="#"/>	<input type="text" value="#"/> <input type="text" value="#"/>	<input type="text" value="#"/> <input type="text" value="#"/>	<input type="text" value="None"/> <input type="text" value="#"/> <input type="text" value="#"/>
<input type="checkbox"/> 3:	<input type="text" value="#"/>	<input type="text" value="#"/> <input type="text" value="#"/>	<input type="text" value="#"/> <input type="text" value="#"/>	<input type="text" value="None"/> <input type="text" value="#"/> <input type="text" value="#"/>
<input type="checkbox"/> 4:	<input type="text" value="#"/>	<input type="text" value="#"/> <input type="text" value="#"/>	<input type="text" value="#"/> <input type="text" value="#"/>	<input type="text" value="None"/> <input type="text" value="#"/> <input type="text" value="#"/>
<input type="checkbox"/> 5:	<input type="text" value="#"/>	<input type="text" value="#"/> <input type="text" value="#"/>	<input type="text" value="#"/> <input type="text" value="#"/>	<input type="text" value="None"/> <input type="text" value="#"/> <input type="text" value="#"/>
<input type="checkbox"/> 6:	<input type="text" value="#"/>	<input type="text" value="#"/> <input type="text" value="#"/>	<input type="text" value="#"/> <input type="text" value="#"/>	<input type="text" value="None"/> <input type="text" value="#"/> <input type="text" value="#"/>
7:	Standard demographics			<input type="text" value="#"/> <input type="text" value="#"/>
8:	Not within expected values			<input type="text" value="#"/> <input type="text" value="#"/>
Critical Limits	Low <input type="text" value="#"/>	High <input type="text" value="#"/>	Unit <input type="text" value="mg/L"/>	Select <input type="text"/> Decimal Places <input type="text" value="2"/>

Calibrators	General	ISE		
Test Name: CYS1G <input type="button" value="Use Serum Cal."/>				
Calibration Type:	<input type="text" value="6AB"/>	Formula: <input type="text" value="Spline"/> Counts: <input type="text" value="2"/>		
<Calibrator Parameters>		Slope Check <input type="text" value="+"/>		
	Calibrator	OD		
		Conc		
		Range		
		Low		
		High		
Point-1	CYSC Calibrator Level 1	†	-2.0000	2.0000
Point-2	CYSC Calibrator Level 2	†	-2.0000	2.0000
Point-3	CYSC Calibrator Level 3	†	-2.0000	2.0000
Point-4	CYSC Calibrator Level 4	†	-2.0000	2.0000
Point-5	CYSC Calibrator Level 5	†	-2.0000	2.0000
Point-6	CYSC Calibrator Level 6	†	-2.0000	2.0000
Point-7				
MB Type Factor	<input type="text"/>	1-Point Calibration Point <input type="text" value="None"/>	<input type="checkbox"/> with Conc-0	Interval (RB) <input type="text"/>
			Stability	Interval (ACAL) <input type="text"/>
			Reagent Blank	<input type="text" value="28"/> Day <input type="text" value="0"/> Hour
			Calibration	<input type="text" value="28"/> Day <input type="text" value="0"/> Hour

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