# Sistatin C

## Beckman Coulter® AU Sistemlerinde (AU5800, AU680, AU480, AU2700, DxC 700 AU) Sistatin C İmmünoassay İçin Kullanma Talimatı

REF B08179

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Sistem Reaktifi: B08179 Reaktif ID: 228

Specific Test Parameters										
General	LIH	ISE	Range							
Test Name:	CysC	<	>	Type:	Serum	Operation:	Yes			
Sample:	Volume	2	µL	Dilution	0	µL	Pre-Dilution Rate:	1		
Reagents:	R1 Volume	150	µL	Dilution	0	µL	Min OD		Max OD	
	R2 Volume	30	µL	Dilution	10	µL	L		H	
Wavelength:	Pri.	540	∇	Sec.		∇	Reagent OD limit:			
Method:	End Point	∇	First L	-2.0	First H	2.0	Last L		Last H	
Reaction slope:	+	∇	Dynamic Range:	L	0.4	H	7.8			
Measuring Point 1:	First	13	Last	27	Correlation Factor:	A	1.00	B	0.00	
Measuring Point 2:	First		Last		On-board stability period:		60			
Linearity:		%								
No Lag Time:		∇								

Specific Test Parameters							
General	LIH	ISE	Range				
Test Name:	CysC	<	>	Type:	Serum		
Value/Flag:	#	∇	Level L:	#	Level H:	#	
Normal Ranges:	Age L	Age H					
	Sex	Year	Month	Year	Month	L	H
<input type="checkbox"/>	1.	#	#	#	#	#	#
<input type="checkbox"/>	2.	#	#	#	#	#	#
<input type="checkbox"/>	3.	#	#	#	#	#	#
<input type="checkbox"/>	4.	#	#	#	#	#	#
<input type="checkbox"/>	5.	#	#	#	#	#	#
<input type="checkbox"/>	6.	#	#	#	#	#	#
<input type="checkbox"/>	7.	None Selected				#	#
<input type="checkbox"/>	8.	Out of Range	L	H		#	#
Panic Value:	#	#	Unit:	mg/L	Decimal places:	#	

Calibration Specific							
General	ISE						
Test Name:	CysC	<	>	Type:	Serum		
Calibration Type:	6AB	Formula:	Spline	Counts:	#	Process:	CONC
Point 1:	Cal. No.	OD	CONC	Factor/OD-L	Factor/OD-H		
Point 2:	1		**				
Point 3:	2		**				
Point 4:	3		**				
Point 5:	4		**				
Point 6:	5		**				
Point 7:	6		**				
1-Point Cal. Point:	<input type="checkbox"/>	with CONC-0	Slope Check:	+	Advanced Calibration:	#	∇
MB Type Factor:		Calibration Stability Period:	28				

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Sistem Reaktifi: B08179

Reaktif ID: 228

Parameters		Specific Test Parameters		
General	LIH	ISE	Calculated Test	Range
Test Name: <input type="text" value="CysC"/> < > Type: <input type="text" value="Serum"/> Operation: <input type="text" value="Yes"/>				
Sample Volume	<input type="text" value="2"/> µL	Dilution	<input type="text" value="0"/> µL	OD Limit
Pre-Dilution Rate	<input type="text" value="1"/>			Min.OD <input type="text"/> Max.OD <input type="text"/>
Rgt. Volume	R1(R1-1) <input type="text" value="150"/> µL	Dilution	<input type="text" value="0"/> µL	Reagent OD Limit
				First Low <input type="text" value="-2.0"/> High <input type="text" value="2.0"/>
				Last Low <input type="text"/> High <input type="text"/>
R2(R2-1)	<input type="text" value="30"/> µL	Dilution	<input type="text" value="10"/> µL	
Common Rgt. Type	<input type="text"/>	Name	<input type="text"/>	Dynamic Range Low <input type="text" value="0.4"/> High <input type="text" value="7.8"/>
Wavelength	Pri <input type="text" value="540"/> nm	Sec. <input type="text"/>		Correlation Factor A <input type="text" value="1.00"/> B <input type="text" value="0.00"/>
Method	<input type="text" value="End Point"/>			Factor for Maker A <input type="text"/>
Reaction Slope	<input type="text" value="+"/> ▾			Onboard Stability Period <input type="text" value="60"/> Day <input type="text"/> Hour
Measuring Point1 First	<input type="text" value="13"/>	Last <input type="text" value="27"/>		LIH Influence Check <input type="text"/>
Measuring Point2 First	<input type="text"/>	Last <input type="text"/>		Lipemia <input type="text"/>
Linearity Limit	<input type="text"/>			Icterus <input type="text"/>
Lag Time Check	<input type="text"/>			Hemolysis <input type="text"/>

Parameters		Specific Test Parameters						
General	LIH	ISE	Calculated Test	Range				
Test Name: <input type="text" value="CysC"/> < > Type: <input type="text" value="Serum"/>								
Value/Flag:	<input type="text" value="#"/>	Low <input type="text"/>	High <input type="text"/>	Panic Value				
Level				Low <input type="text"/> High <input type="text"/>				
Specific Ranges:	From	To						
<input type="checkbox"/> 1.	Sex <input type="text"/>	Year <input type="text"/>	Month <input type="text"/>	Year <input type="text"/>	Month <input type="text"/>	Low <input type="text"/>	High <input type="text"/>	
<input type="checkbox"/> 2.	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	
<input type="checkbox"/> 3.	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	
<input type="checkbox"/> 4.	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	
<input type="checkbox"/> 5.	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	
<input type="checkbox"/> 6.	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	
<input type="checkbox"/> 7.	No demographics						<input type="text"/>	<input type="text"/>
<input type="checkbox"/> 8.	Not within expected values						<input type="text"/>	<input type="text"/>
Unit	<input type="text" value="mg/L"/>	Decimal Places	<input type="text" value="#"/>					

Parameters		Calibration Parameters		
Calibrators	Calibration Specific	STAT Table Calibration		
General	ISE			
Test Name: <input type="text" value="CysC"/> < > Type: <input type="text" value="Serum"/> <input type="checkbox"/> Use Serum Cal.				
Calibration Type:	<input type="text" value="6AB"/>	Formula:	<input type="text" value="Spline"/>	Counts: <input type="text" value="#"/>
<Calibrator Parameters>				
Calibrator	OD	Conc	Factor Range	Slope Check <input type="text" value="+"/>
Point 1:	<input type="text" value="1"/>	<input type="text"/>	<input type="text"/>	
Point 2:	<input type="text" value="2"/>	<input type="text"/>	<input type="text"/>	Allowance Range Check
Point 3:	<input type="text" value="3"/>	<input type="text"/>	<input type="text"/>	
Point 4:	<input type="text" value="4"/>	<input type="text"/>	<input type="text"/>	<input type="checkbox"/> Reagent Blank <input type="text"/>
Point 5:	<input type="text" value="5"/>	<input type="text"/>	<input type="text"/>	<input type="checkbox"/> Calibration <input type="text"/>
Point 6:	<input type="text" value="6"/>	<input type="text"/>	<input type="text"/>	Advanced Calibration
Point 7:	<input type="text"/>	<input type="text"/>	<input type="text"/>	Operation <input type="text" value="#"/>
Point 8:	<input type="text"/>	<input type="text"/>	<input type="text"/>	Interval (RB/ACAL) <input type="text" value="#"/>
Point 9:	<input type="text"/>	<input type="text"/>	<input type="text"/>	
Point 10:	<input type="text"/>	<input type="text"/>	<input type="text"/>	
<Point Cal. For Master Curve>	No. of Correction Points <input type="text"/>	Use Master Curve <input type="checkbox"/>	Lot Calibration <input type="checkbox"/>	
Calibrator	OD	Conc	Factor Range	Stability
Point-1	<input type="text"/>	<input type="text"/>	<input type="text"/>	Reagent Blank <input type="text" value="28"/> Day <input type="text"/>
Point-2	<input type="text"/>	<input type="text"/>	<input type="text"/>	Calibration <input type="text" value="28"/> Day <input type="text"/>
MB Type Factor:	<input type="text"/>	1-Point Calibration Point <input type="checkbox"/>	<input type="text"/>	<input type="checkbox"/> with Conc-0

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Sistem Reaktifi: B08179  
Test adı: CYS1G

Reaktif ID: 228

General	LIH	ISE	Calculated Test	Range
<b>Test Name:</b> CYS1G <b>Test No</b> <b>Type:</b> Serum <b>Operation</b> Yes				
Sample Volume	2.0 µL	Dilution	0 µL	OD Limit
Pre-Dilution Rate	1			Min. OD <input type="text"/> Max OD <input type="text"/>
Reagent Volume	R1 (R1-1) 150 µL	Dilution	0 µL	Reagent OD Limit 1 <sup>st</sup> Low -2.0000 High 2.0000
	R1-2 <input type="text"/> µL	Dilution	<input type="text"/> µL	Last Low -2.0000 High 2.0000
	R2 (R2-1) 30 µL	Dilution	10 µL	Analytical Measuring Range Low 0.4 High 7.8
Common Reagent	Type None	Name	None	Correlation Factor A 1 B 0
Wavelength	Pri 540 nm	Sec	None nm	Manufacturer Factor A 1 B 0
Method	END			
Reaction Slope	+			Onboard Stability Period 60 Day 0 Hour
Measuring Point-1	1st 13	Last	27	LIH Influence Check No
Measuring Point-2	1st <input type="text"/>	Last	<input type="text"/>	Lipemia <input type="text"/>
Linearity Limit	<input type="text"/> %			Icterus <input type="text"/>
Lag Time Check	<input type="text"/>			Hemolysis <input type="text"/>

General	LIH	ISE	Calculated Test	Range
<b>Test Name:</b> CYS1G <b>Test No</b> <b>Type:</b> Serum				
Value/Flag	Value	Level	Low -99999.99	High 99999.99
<b>Specific Ranges</b>				
	Sex	Year	Month	Year
From				
To				
Other Type				
Low				
High				
7: Standard demographics				
8: Not within expected values				
Critical Limits	Low #	High #	Unit mg/L	Select
				Decimal Places 2

Calibrators	General	ISE
<b>Test Name:</b> CYS1G <b>Type:</b> Serum		
<input type="checkbox"/> Use Serum Cal.		
Calibration Type:	6AB	Formula: Spline
		Counts: 2
<Calibrator Parameters>		
		Slope Check +
		Allowable Range Check
		<input type="checkbox"/> Reagent Blank
		<input type="checkbox"/> Calibration
		Advanced Calibration
		Operation No
		Interval (RB)
		Interval (ACAL)
		Stability
		Reagent Blank 28 Day 0 Hour
		Calibration 28 Day 0 Hour
MB Type Factor		1-Point Calibration Point None
		<input type="checkbox"/> with Conc=0

	Calibrator	OD	Conc	Range	
				Low	High
Point-1	CYSC Calibrator Level 1		†	-2.0000	2.0000
Point-2	CYSC Calibrator Level 2		†	-2.0000	2.0000
Point-3	CYSC Calibrator Level 3		†	-2.0000	2.0000
Point-4	CYSC Calibrator Level 4		†	-2.0000	2.0000
Point-5	CYSC Calibrator Level 5		†	-2.0000	2.0000
Point-6	CYSC Calibrator Level 6		†	-2.0000	2.0000
Point-7					

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